

NATURA[®]

Accordion Flange



Flat Base Plates with Accordion Flange[®] and Convatec Mouldable Technology[™]

Stomahesive[®] Accordion Flange[®] base plate with hydrocolloid collar



Product Code		Base Plate Size	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
421033	●	57mm	13-22mm	Mouldable	10	80025W	30
421034	●	57mm	22-33mm	Mouldable	10	80025W	30
421035	●	70mm	33-45mm	Mouldable	10	80025W	30

Durahesive[®] Accordion Flange[®] base plate with hydrocolloid collar



Product Code		Base Plate Size	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
421039	●	57mm	13-22mm	Mouldable	10	80025W	30
421040	●	57mm	22-33mm	Mouldable	10	80025W	30
421041	●	70mm	33-45mm	Mouldable	10	80025W	30

FEATURES

- The accordion flange lifts to provide generous finger room so coupling is easy and comfortable.
- Accordion Flange® has a low profile in the 'DOWN' position so it's not visible under clothing. The base plate moves as the body flexes
- Available in flat cut to fit, Convex cut to fit and flat ConvaTec Mouldable Technology™
- Connects with an audible 'click' and physical 'snap' to confirm the pouch is firmly attached

Designed for patients

- with dexterity issues as it provides finger room so coupling is easy and comfortable
- to prevent pressure or trauma if abdominal tenderness is present due to post-operative stoma surgery/ incisions

HOW IT WORKS

